



The safety of the reciprocation movement reaching a higher level of efficiency

M-ONE GOLD		21mm	25mm	31mm
Shaping File				
	20 .07	MOG2120	MOG2520	MOG3120
	25 .07	MOG2125	MOG2525	MOG3125
	35 .06	MOG2135	MOG2535	MOG3135
	45 .05	MOG2145	MOG2545	MOG3145
Glidepath File				
	15 .02	MOG2115	MOG2515	MOG3115



Sold In Packs of 3
\$20.00




For Information Contact Us
Ph. 800.221.0750
Or Visit: www.Medidenta.com

Suggested use (based on average morphology):

Curved and narrow canals

Lower molars MB, MB2 and ML
Upper molars MB and DB



CS1 18 .02


CS2 20 .04

AF1 20 .07

Optional M-File Gold
35 .04 in rotation

Straight canals

Lower premolars
Anterior



20 .07


35 .06

45 .05

Optional M-File Gold
larger files in rotation

Curved canals

Upper premolars
Lower molar D
Upper molar P



AF2 25 .08

AF3 30 .09

Optional M-File Gold
40 .04 in rotation